



Genemed Synthesis, Inc.

6203 Wood Lake Center Dr., Bldg. 2, San Antonio, TX 78244, USA
Toll free (800) 344-5337; **Phone:** (210) 745-5988; **Fax** (210) 745-5992
Email: info@genemedsyn.com, **Website:** www.genemedsyn.com

Super Fluor™ 555 acid

Cat #	FD-1006	Unit Size	5 mg
MW			
Ex (nm)	556	Em (nm)	567
Solvent	DMF	Storage	F/D/L

Cy3 is the preferred dye for preparing orange fluorescent bioconjugates. Compared to the spectra of Cy3 dyes, the spectra of Superfluor 555 conjugates are slightly red-shifted, resulting in an optimal match to filters designed for Cy3 dyes. Superfluor 555 is more photostable than Cy3, providing researchers with additional time for image capture. Additionally, our in-house data indicated that Superfluor 555 is superior to Cy3 in fluorescence polarization-based assays.

Rev. 80725F