



Genemed Synthesis, Inc.

6203 Wood Lake Center Dr., Bldg. 2, San Antonio, TX 78244, USA
Toll free (800) 344-5337; **Phone:** (210) 745-5988; **Fax** (210) 745-5992
Email: info@genemedsyn.com, **Website:** www.genemedsyn.com

Super Fluor™ 680 succinimidyl ester

Cat #	FD-1035	Unit Size	1 mg
MW			
Ex (nm)	680	Em (nm)	702
Solvent	DMF	Storage	F/D/L

With a peak excitation at 678 nm and emission maximum at 699 nm, Superfluor 680 dyes are spectrally similar to Cy5.5 dyes. Fluorescence emission of Superfluor 680 dyes is well separated from that of other commonly used red fluorophores, such as TAMRA, R-phycoerythrin and Superfluor 647 dyes, making it ideal for three- and four-color labeling.

Rev. 80725F